



SparSan Q™

disinfectant deodorant

SparSan Q Disinfectant Deodorant is a versatile, ready to use aerosol disinfectant deodorant. It is ideal for convenient disinfection after spot clean up, and for those areas not easily reached during routine disinfection procedures. Because it's an aerosol, there is no mixing or any waste.



Providing Quality Maintenance and Chemical Specialty Solutions with Unparalleled Customer Service



Spartan Chemical Company, Inc.
1110 Spartan Drive
Maumee, OH 43537
1-800-537-8990
www.spartanchemical.com

- Kills HIV-1 on Pre-Cleaned Environmental Surfaces
- Eliminates odors and has a linen clean fragrance
- Germicidal
- Fungicidal
- Pseudomonacidal
- Tuberculocidal
- Kills Household Germs
- Prevents Mold and Mildew

HOSPITAL GRADE DISINFECTANT:

Disinfects and deodorizes surfaces through effective bactericidal and fungicidal activity. Kills *Staphylococcus aureus* (Staph), *Streptococcus hemolyticus* (Strep), *Pseudomonas aeruginosa*, *Mycobacterium bovis* (BCG), *Salmonella choleraesuis* and *Trichophyton mentagrophytes* (athletes foot fungus) on hard, non-porous environmental surfaces. Meets EPA Efficacy Standards for hospital aerosol disinfectants.

SPRAYED ON SURFACES:

Deodorizes by disinfecting. Germicidal activity kills most odor-producing organisms. Controls unpleasant odors. Pleasantly scented. For use in hospital rooms, sick rooms, doctors' examining rooms, kitchens, nursery and bathrooms.

EXCELLENT DEODORIZER:

SparSan Q Disinfectant Deodorant destroys bacteria the initial source of odors. It is not a cover up. Deodorizers which do not have disinfectant properties merely cover one odor with another.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. SHAKE WELL BEFORE USING. Prior to application, thoroughly clean all surfaces to be treated to remove gross filth or heavy soil.

To Disinfect: Hold can upright and spray from a distance of 6 to 8 inches until area for treatment is thoroughly moistened and allow to dry for 10 minutes before wiping or allow to air dry.

For Control of Mold and Mildew on Hard, Non-Porous Environmental Surfaces: Spray surface until thoroughly moistened. Allow to air dry. Repeat application in 7 days.

For Control of Fungus in areas such as bathrooms, locker rooms, floors, benches, shower stalls or other non-porous environmental surfaces commonly contacted by bare feet: Spray surface until thoroughly moistened. Allow to air dry.

For Use On: Walls, floors, bed frames, bedpans, urinals, fixtures, containers, toilet seats, diaper pails, telephones, trash cans, and sickrooms. Due to the possibility of electric shock when used on telephones, avoid excess spraying of telephone.

Kills Hiv On Pre-Cleaned Environmental Surfaces And Objects Previously Soiled With Blood/Body Fluids in health care settings. Or in other settings where there is a likelihood of soiling of inanimate surfaces or objects with blood or body fluids in which the surfaces and objects can be associated with the potential for transmission of human immunodeficiency virus type 1 (HIV-1) (associated with AIDS).

Special instructions for cleaning and decontamination against hiv on surfaces and objects soiled with blood or body fluids:

Personal protection:

Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks or eye coverings.

Cleaning Procedure: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of disinfectant.

Disposal Of Infectious Material: Blood and other body fluids should be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal.

Contact Time: Leave surfaces wet for 10 minutes. This contact time will not control other common types of viruses and bacteria.

Be sure to read all Directions, Precautionary and First Aid Statements on product labels before use of this or any Spartan product. If questions remain, consult your employer or a physician. Material Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor or by visiting www.spartanchemical.com.

SPECIFICATION DATA:

pH – 10.5 – 11.5

Flash Point - Flammable

Stability - Shelf @ 24°C/75°F - One year minimum



ISSA MEMBER
The Experts
on Cleaning and Maintenance

Distributed by:

PACKAGING

607600 20-ounce cans (12 per case)



Label copy is provided in English and Spanish.

GUARANTEE: Spartan's modern manufacturing and laboratory control insure uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within one year of the date of manufacture. Use product as directed and read all precautionary statements.

For institutional and industrial use only.

© SCC 5/11 L6076